

Wash Baskets & Machine Accessories

For use in Miele Large Chamber Washers Models PG 8527 & PG 8528





Miele Baskets & Accessories: Advancing the Science of Clean

Engineered Simplicity

Miele washers expertly clean a wide variety of items through a simple-to-use system of interchangeable baskets and inserts. Miele wash baskets are engineered to roll directly onto the transfer cart, with a modular basket design allowing simple removal of each level away from the machine. This ergonomic design allows easy loading and unloading at the washer or at the glassware storage cabinet.

This interchangeable basket system ensures a precise match to your specific cleaning requirements. And, thanks to Miele's universal design for large chamber wash baskets and inserts, a new basket will fit a new Miele PG 8527 washer as precisely as it will an older model.

To further simplify basket and wash cycle selection, the PG 8527 features an optional Automatic Basket Recognition, ensuring repeatable wash results and eliminating guesswork from the machine operator regarding wash parameters.

Injection Cleaning and Drying

The Miele PG 8527 washer features 1 to 3 levels of injection cleaning or up to 5 levels of non-injection cleaning, providing excellent throughput and wide array of wash load possibilities.

Miele wash baskets feature removable and interchangeable injectors for a truly custom basket configuration. This interchangeability also allows easy replacement of older injectors, ensuring many years of useful operation.

When Miele's HEPA-filtered forced air drying is selected, air is circulated throughout the wash chamber and through the injectors to the interior of narrow-neck flasks. This provides superior drying typically within 15 to 30 minutes, eliminating the need for a separate drying oven. This system frees up valuable floor space, reduces labor costs, and lowers your electrical costs.

Gentle Yet Effective Cleaning

Miele washers utilize high-volume, low-pressure water circulation to prevent breakage inside the wash chamber.

This care is carried over into our wash basket design. Each basket delicately cradles the item being washed while providing renowned, repeatable wash results. An excellent example is pipettes. Miele provides perfect results for difficult items like this, eliminating breakage even for small items.

Solutions for Your Specific Application

If you can not find a basket solution, please contact us for more information. Through the exclusive Miele Application Laboratory, housed in Princeton, NJ, with a sister lab in Gütersloh, Germany, we are continually developing new wash baskets, ensuring you a competent, customized solution to your cleaning requirements.

Mobile Baskets



Spray Arm Monitoring System

These Mobile Units include spray arm monitoring for each spray arm on the basket. This system ensures that all spray arms rotate at the proper speed and creates an alarm condition if a spray arm is blocked due to improper loading.

This is another example of how Miele ensures valid wash results.



E 975/1 Mobile Unit

- 2 levels w/ built-in spray arm
- Magnetic coding strip for automatic sensing of a specific wash routine
- Clearance: (Level 1 = Lower Level)
 - 1: 297 mm H, 592 mm W, 780 mm D
 - 2: 290 mm H, 592 mm W, 780 mm D
- Dimensions: Height 427 mm
Width 640 mm
Depth 790 mm
- Article No. 69597502



E 935/1 Mobile Unit

- 3 levels w/ 2 built-in spray arms
- Magnetic coding strip for automatic sensing of a specific wash routine
- Clearance: (Level 1 = Lower Level)
 - 1: 202 mm H, 590 mm W, 780 mm D
 - 2: 202 mm H, 590 mm W, 780 mm D
 - 3: 132 mm H, 590 mm W, 780 mm D
- Dim.: 524 mm H, 640 mm W, 790 mm D
- Article No. 69593502



E 900-4/1 Mobile Unit

- 4 levels w/ 3 built-in spray arms
- Magnetic coding strip for automatic sensing of a specific wash routine
- Clearance: (Level 1 = Lower Level)
 - 1: 112.5 mm H, 593 mm W, 780 mm D
 - 2: 112.5 mm H, 593 mm W, 780 mm D
 - 3: 112.5 mm H, 593 mm W, 780 mm D
 - 4: 114 mm H, 593 mm W, 780 mm D
- Dimensions: Height 557 mm
Width 640 mm
Depth 790 mm
- Article No. 69590003



E 900-5/1 Mobile Unit

- 5 levels w/ 4 built-in spray arms
- Magnetic coding strip for automatic sensing of a specific wash routine
- Clearance:
 - 1: 80 mm H, 593 mm W, 780 mm D
 - 2: 80 mm H, 593 mm W, 780 mm D
 - 3: 80 mm H, 593 mm W, 780 mm D
 - 4: 80 mm H, 593 mm W, 780 mm D
 - 5: 73 mm H, 593 mm W, 780 mm D
- Dim.: 605 mm H, 640 mm W, 790 mm D
- Article No. 69590004

Direct Injection Baskets



E 940 Mobile Injector

- With two levels for direct injection washing and drying of narrow-necked laboratory glassware
- Upper level can be removed for single level cleaning of tall items
- Magnetic coding strip for automatic sensing of a specific wash routine
- Connection for forced air drying through injectors
- 115 injectors with clips are removable for custom configuration of the basket
- Lower level = 35 injectors (E 352), 35 clips (E 354), 6 x 220 mm
- Upper level = 80 injectors (E 351) 80 clips (E 353), 4 x 160 mm
- Dimensions: Height 565 mm
Width 640 mm
Depth 790 mm
- Article No. 69594001



E 957 Mobile Injector Unit

- For 1-12 large items (carboys)
- Includes twelve injector jets and twelve plastic coated supports to secure glassware
- Height-adjustable holding frame with 3 short and 2 long struts for various glass diameters
- Magnetic coding strip for automatic sensing of a specific wash routine
- Dimensions: Height 353 mm
Width 640 mm
Depth 790 mm
- Article No. 69595701



E 950/1 Mobile Injector

- With three levels for direct injection washing and drying of narrow-necked laboratory glassware
- Middle level can be removed for cleaning of tall items
- Magnetic coding strip for automatic sensing of a specific wash routine
- Connection for forced air drying through injectors
- 232 injectors with clips are removable for custom configuration of the basket
- Levels 1 & 3 = 80 injectors
- Level 2 = 72 injector
- Clearance: (Level 1 = Lower Level) 1/3: 148 mm H
- Dimensions: Height 572 mm
Width 640 mm
Depth 790 mm
- Article No. 69595002



Modular Baskets: Design Your Own
Miele's modular basket concept allows the user to design a unique wash basket through the combination of various inserts, connected by a simple docking system.

When configured, the wash basket can be a combination of only injection cleaning, only non-injection cleaning, and any variation in-between.

This system allows single level cleaning of tall items, dual level cleaning of smaller items, or split-level cleaning, allowing both tall items and small items can be cleaned in a single wash load.

For even greater throughput, the modular basket configuration provides an easy load, unload, and transfer. Unlike other manufacturers, the G 7827 modular basket allows the operator to remove one insert at a time, eliminating the need to load and unload at the machine or lift the entire large, heavy basket.

Modular Baskets: Design Your Own



E 941 Modular Mobile Unit

- For direct injection cleaning, standard cleaning, or a combination of both in a single load on one or two levels
- 2 modular inserts possible each level
- Magnetic coding strip for automatic sensing of a specific wash routine
- Forced-air drying through injectors
- Dimensions: Height 421 mm
Width 619 mm
Depth 790 mm
- Article No. 69594101



E 947 Modular Injector Insert

- For centrifuge tubes, vials & test tubes
- For use with E 941
- Includes 88 jets, 110 mm L / Ø 2.5 mm
- Dimensions: Height 170 mm
Width 558 mm
Depth 352 mm
- Article No. 69594701



E 942/1 Modular Injector Insert

- Injector unit for 116 pipettes
- For use with E 941
- Forced-air drying through injectors
- Each section is 16 x 16 mm
- Frame supports wide mouth glassware
- Dimensions: Height 279 mm
Width 558 mm
Depth 352 mm
- Article No. 69594202



E 945 Modular Insert Support Frame

- Support frame for use with E 941 for holding inserts for wide mouth glassware such as the E 960, E 963, and E 965
- Dimensions: Height 55 mm
Width 558 mm
Depth 352 mm
- Article No. 69594501



E 943 Modular Injector Insert

- Injector unit for narrow-necked laboratory glassware, 100-500 ml
- For use with E 941
- Forced-air drying through injectors
- Injectors are removable for custom configuration
- 32 Injectors (E 351), 4 x 160 mm
- 32 Mielava injector clips (E 353)
- Dimensions: Height 190 mm
Width 558 mm
Depth 352 mm
- Article No. 69594301



E 944 Modular Injector Insert

- Injector unit for narrow-necked laboratory glassware, 500-1000 ml
- For use with E 941
- Forced-air drying through injectors
- Injectors are removable for custom configuration
- 15 injectors (E 352), 6 x 220 mm
- 15 Mielava injector clips (E 354)
- Dimensions: Height 250 mm
Width 558 mm
Depth 352 mm
- Article No. 69594401

Modular Non-Injection Baskets



- E 945 Modular Insert Support Frame**
- Support frame for use with E 941 for holding inserts for wide mouth glassware such as the E 960, E 963, and E 965
 - Article No. 69594501



- E 984 Stainless Steel Basket**
- For various items
 - E 986 Stainless Steel Support available for tall objects such as trays
 - Dimensions: Height 65 mm
w/ handles 105 mm
Width 371 mm
Depth 572 mm
 - Article No. 69598401

- E 986 Stainless Steel Support**
- Article No. 69598601



- Stainless Steel Basket with Lid**
- Ideal for larger loose items
 - Article No. 78270034



- E 969 Mesh Tray**
- For various items
 - 2 A 19 Lids needed for complete coverage of E 969 (Lids not included)
 - Glassware not included
 - Dimensions: Height 67 mm
w/ handles 122 mm
Width 363 mm
Depth 533 mm
 - Article No. 69596901

- A 19 Lid (for E 969)**
- Article No. 69796901



E 960/1 Insert

- Contains 20 large and 26 small spring clips to secure glassware
- For wide-necked Erlenmeyer flasks and measuring flasks
- Dimensions: Height 185 mm
Width 357 mm
Depth 522 mm
- Article No. 69596002



E 963 Insert

- Contains 33 holders
- For beakers up to 250 ml or other wide mouthed items not requiring injection
- For use in E 935, E 975, or E 945
- Dimensions: Height 155 mm
Width 357 mm
Depth 522 mm
- Article No. 69596301



E 965 Insert

- Contains 15 holders
- For beakers 250-600 ml or other wide mouthed items not requiring injection
- For use in E 935, E 975, or E 945
- Dimensions: Height 173 mm
Width 357 mm
Depth 522 mm
- Article No. 69596501



Plastic-Coated Rack for Glass Plates

- Coating provides greater protection to glass plates
- Article No. 78270031

Machine Accessories: PG 8527, PG 8528, G 7827 & G 7828



SBW Base

- For PG 8527 & PG 8528
- Stationary base
- Article No. 69253001

SBW/3 Base

- For G 7827 & G 7828
- Stationary base
- Article No. 69371003



TAE HEPA-filtered Drying Unit

- For PG 8527 & PG 8528
- Article No. 69250001

TAE HEPA-filtered Drying Unit

- For G 7827 & G 7828
- Article No. 69201001



MAV 27/28 Steam Condenser Enclosure

- For PG 8527 & PG 8528
- Required with DK 27/28
- Article No. 69252001

MAV 27/28 Steam Condenser Enclosure

- For G 7827 & G 7828
- Required with DK 27/28
- Article No. 69210002



SBWR Roller Base

- For PG 8527 & PG 8528
- Recommended if machine must be installed <18" from wall
- Simplifies machine service
- Article No. 69371005

SBW/R Roller Base

- For G 7827 & G 7828
- Recommended if machine must be installed <18" from wall
- Simplifies machine service
- Article No. 69371005



DK 27/28 Steam Condenser

- For PG 8527 & PG 8528
- Requires MAV 27/28
- Article No. 69251001

DK 27/28 Steam Condenser

- For G 7827 & G 7828
- Requires MAV 27/28
- Article No. 69230004



MF 27/28 Transfer Cart

- For PG 8527 & PG 8528
- For safe and easy transfer of baskets
- Docks & couples with washer
- Article No. 69200110

MF-1 Transfer Cart

- For G 7827 & G 7828
- For safe and easy transfer of baskets
- Docks & couples with washer
- Article No. 69200106



Indicator Light Box

- Indicates progress of wash cycle

3 Light Status Indicator Box

- Shown above
- Article No. PROF18975T999

5 Light Status Indicator Box

- Article No. 78270022



PRT/1 Printer

- Shown above
- Article No. 69247001

External Printer

- Article No. 5035-033-01

PG 8527/PG 8528 Trim Kit

- For in-wall installations
- 3" wide stainless steel trim to cover wall opening most commonly for pass-through installations
- Article No. PROFTRIMKITPG8527

G 7827/G 7828 Trim Kit

- For in-wall installations
- Article No. PROFTRIMKITG7827





Machine Accessories: G 7895/1 Pure Water System



G 7895/1 Pure Water System

- Stainless steel (shown) or white cabinet
- Inc. conductivity meter (in microS/cm)
- Inc. (1) E 318 cartridge & (1) E 315 resin
- Unit holds max. of two E 318 cartridges
- Dimensions: Height 850 mm
w/out lid 820 mm
Width 300 mm
Depth 600 mm
- Two pressure hoses, approx. 4' long
- Water intake to 3/4" Male Hose Thread
- Min. 36 psi incoming water pressure

G 7895/1 - Stainless Steel

- Article No. 69789504

G 7895/1 - White

- Article No. 69789503



E 318 Stainless Steel Cartridge

- Stainless steel cartridge
- Pressure-proof w/ vent & relief valve
- Holds 20 Liter of E 315 disposable resin
- Diameter 240 mm, Height 570 mm
- Article No. 69431801



E 316 Refill Kit for Resin Replacement

- Refill kit for resin replacement
- Inc. 30 Liter plastic barrel and funnel
- Capacity: 12 gallon/45 Liter
- Article No. 69431601



E 315 Disposable Resin

- Disposable resin, 20 Liter
Two 10 Liter bags
- Vacuum-packed plastic bags
- Inc. filter bag for replacement
- Article No. 69431501



Miele
PROFESSIONAL

Miele
PROFESSIONAL



Service and Validation



Validation Documentation

- Complete IQ/OQ
- User-friendly, step-by-step check-offs
- Covers all alarm conditions
- Assures safe validation
- Assures documents are up-to-date



Validation Execution

- Miele-direct technicians
- Very professional and expertly trained
- Prompt
- Set, competitive pricing alleviates cost overruns



Preventive Maintenance Plans

- Assures you of maximum machine performance and longevity
- Various plans to meet your requirements and budgets
- Performed by Miele direct or factory trained technicians

Miele Professional

9 Independence Way
Princeton, NJ 08540
USA

& 800.991.9380

☎ 609.419.4241

🌐 labwasher.com

Miele Limited

161 Four Valley Drive
Vaughan, ON L4K 4V8
Canada

& 800.643.5381

☎ 905.532.2290

🌐 miele.ca

Miele S.A. de C.V.

Av. Santa Fe 170, German Centre: 0-4-2
Col. Lomas de Santa Fe
C.p. 01210. Mexico, D.F.

& (+5255) - 8503 9870 al 73

☎ (+5255) - 8503 9874

🌐 miele-professional.com.mx

Every effort has been made to ensure that the information given in this brochure is correct at the time of printing. However, due to continuous product improvement, Miele reserves the right to make changes to products and technical data without prior notice.