

EA-34004C SMARTBOX™ Lab Control Unit

Operates at a fixed flow rate based on the external CO₂ source.

The Controller automates on/off gas flow at the push of a button, delivering CO₂ to individual cages or Euthanex lids. The EA-34004C has four individual gas outputs allowing up to four cages to be euthanized individually or staggered. Increasing the user efficiency when large of number of animals need to be euthanized individually. The EA-34004C has four separate switches that offer eight timing presets. Presets are programmable. This economical unit is designed for small scale use by individual investigators in their laboratories.



Lids and chambers not included.

Lab Control Unit Features:

- AVMA Compliant, introduces gas 10%-30% displacement
- Treats animals at individual cage level
- Designed for investigators that need to retrieve animals quickly
- Utilizing the full range of Euthanex Cage Lids or existing induction chamber
- Increase efficiency staggering cage start times
- Minimizes operator error and excessive gas use
- Automatic shut-off saves gas and protects personnel
- Double throw switches allow for two different time settings for each output
- Easily attaches to wall